

Working together to maintain health and support clean sport

Testing Locations



- USADA screens each test assignment location for COVID-19-related restrictions (e.g., shelter-in-place orders).
- We assign DCOs with priority given to those located closest to the athlete.
- If a DCO arrives and you would like to be tested elsewhere, please communicate that preference.

Monitoring DCOs Health




- USADA DCOs and other sample collection personnel are regularly tested for COVID-19.
- All DCOs complete a health monitoring questionnaire the day they are assigned to test that will further determine their clearance to conduct the assignment.

Using Personal Protective Equipment



- All DCOs will wear a mask throughout each test session.
- All DCOs will wear a face shield after notification has occurred.
- All DCOs will wear gloves when appropriate and available.
- You can wear your own mask, or USADA will provide you with one if needed.

Practicing Social Distancing



- All DCOs will maintain a 6-foot distance whenever possible.
- Please help us maintain the 6-foot distance.

Implementing Mandatory Hand Washing



- Athletes must wash hands with soap and water after notification and sample provision.
- All DCOs must wash hands with soap prior to the sample collection process.
- If soap is not available for any reason,
 DCOs carry hand sanitizer as an alternative hand washing solution.

Disinfecting Surfaces and Materials



 DCOs will bring and use disinfectant wipes and/or spray throughout and following the testing session.

Contact USADA directly with questions or concerns at (719) 785-2000 or (866) 601-2632, and at athleteexpress@USADA.org

Whereabouts Reminders

RTP

Remove and add competitions as events are updated



CAP

Remove and add competitions as events are updated

Update training and school locations as they close and re-open



Update training and school locations as they close and re-open

Provide an accurate primary overnight location



Provide an accurate primary overnight location

Designate a 60-minute window that aligns with your testing preferences

